

# Mobile Calibration System

Model : CalPro



## CalPro

The calibration system developed by Rossmax can check the accuracy of the unit's pressure reading, recalibrate the device when the reading is out of calibration and check / reset the usage record of the unit. The CalPro can also serve as a pressure gauge.



### What's the Mobile Calibration System

Most regulatory agencies suggest recalibration of a unit's pressure readings, by due to the technical constraints, the unit can only be recalibrated by manufacturer or via sophisticated equipment.

Knowing the importance of recalibration to the units and the users, Rossmax devoted 2 years to develop the world's first mobile calibration system. Simply connect Rossmax unit to the system, with its intuitive operation, now recalibration is only a matter of seconds.

Model	Qty per carton	Carton volume
CalPro	4 pcs	0.047 cbm/ctn

○ Storage and Transportation Condition  
-10°C~60°C(14°F~140°F)  
10%~90% RH, 700~1060 hPa

www.rossmax.com  
Epaper\_CalPro\_EN\_ver1346