




Piston Nebulizer


Model: NL100



Patent No.: ZL 2020 3 0115746.3
Other patents pending

Nebulizer Kit
Consisting of:
• Bottle (N2)
• Mouthpiece

- MMAD $\leq 3.8 \mu\text{m}$
- Fine Particle Dose (FPD): 70~75%
- Neb Rate $\geq 0.20 \text{ ml/min}$
- Consistent fine particle size for efficient respiratory treatment
- Mouthpiece and masks for adult and child included
- Compact design



Efficient Treatment



Nebulizer kit




Adult and Child Masks



Aerosol performance tested with 2.5% NaF solution by Cascade Impactor	
Cascade Impactor	
MMAD	FPD
3.390 μm	73.21%

FPD: Fine particle dose, the percentage of particle size less than 5.0 μm



Model	Qty per carton	Carton volume
NL100	8 pcs	0.042 cbm/ctn
QTY per 20'	QTY per 40'	QTY per 40'HQ
5336 pcs	10672 pcs	12384 pcs

Storage and Transportation Condition
-20°C to 60°C (-4°F to 140°F)
10% ~ 90% RH, 700~1060hPa